


## GENERAL FEATURES

- "SAN JOR"
- 304 18/8 stainless steel interior
- Enameled steel exterior
- Glass wool insulation
- 1 Lateral door
- 1 Glass inner door
- Door supplied with a silicone frame
- Adjustable fresh air intake at the top
- Operating cut off switch
- Digital control of temperature Auto-tuning PID BLAST System
- Double viewfinder of internal and programmed temperature
- Temperature range: from 0 to 80 °C (ambient + 5 to 80 °C)
- Light indicators of Resistance, Auto tuning and Alarm
- Accuracy: +- 0,1 °C
- Uniformity: +- 0,2% of the final value of temperature
- Connection 220 Vac. 50-60Hz. Optional 110 Vac
- One year warranty
- Free 2-Year Extended Warranty
- Accessories: 2 adjustable shelves grid type (5 positions)



Inside measurements  
**30 x 40 x 30**



Laboratory Incubator

## TECHNICAL SPECIFICATIONS

<b>Interior</b>	
Width x Height x Depth cm	30 x 40 x 30
Shipping volume (dm <sup>3</sup> )	36
<b>Exterior</b>	
Width x Depth x Height of oven	46 x 45 x 70
Width x Depth x Height with packaging	48 x 48 x 73
<b>Others</b>	
Consumption in watts	130
Net weight Kg	19
Shipping weight Kg	21

- ☺ High performance for permanent use
- ☺ Exclusive PID Control BLAST System
- ☺ Easy use of Digital Control
- ☺ Optimized Auto tuning from factory
- ☺ Low consumption
- ☺ High accuracy
- ☺ Double warranty

[www.sanjor.com.ar](http://www.sanjor.com.ar)